

# INSTRUCTIONS FOR MONITORING HEAD AMP DATA

## Setup

1. Remove the existing head amp data connections and connect the supplied cables as shown in Figure 1. Connect the USB connector from the EZ-Tap to a free USB port on the PC.




Figure 1: Cable Connections

2. Install the EZ-View software, *EZView\_1\_6\_9.msi*, supplied on the USB stick or download from <https://www.stratusengineering.com/downloads/>
3. Start the EZ-View application and select *File->Load Settings* from the menu as shown in Figure 2. Browse to file *Head Amp.evs* supplied on the USB stick. This loads the settings for capturing head amp data. Then select *Settings->Display Mode->Line Mode* from the menu as shown in Figure 3.
4. Start capturing the data by clicking on the green play button in the tool bar. The head amp data should now appear in the window as shown in Figure 4.

5. Data can be captured to a file by selecting File->Save Data as Text or File->Save Data as Binary from the menu.
6. Refer to the following section for an explanation of the data structure.
7. Note that when plugged into the PC the EZ-Tap device creates 2 serial ports; this is so that it can monitor both transmitted and received data. The serial port numbers, e.g. COM4, COM5 are not explicitly referred to in the EZ-View software.
8. The user can then write their own software to capture the data on the appropriate port.




Figure 2: EZ-View Loading Settings




Figure 3: Display Mode

Line	MM/DD/YY hh:mm:ss.ms.us	DTE(hex)	DCE(hex)
8	09/23/20 13:50:39.177.800	ff ff c8 8f ff ff c8 9f	
9	09/23/20 13:50:39.182.600	ff ff c8 af ff ff c8 bf	
10	09/23/20 13:50:39.191.000	ff ff c8 cf 00 00 37 50	
11	09/23/20 13:50:39.199.300	00 ff ff ff ff	
12	09/23/20 13:50:39.208.400		35
13	09/23/20 13:50:39.248.100	10 1b 09 cc 20 1e 74 74	
14	09/23/20 13:50:39.256.500	30 1e 01 8c 40 13 68 1c	
15	09/23/20 13:50:39.273.400	5f 7f ff c8 6f ff ff c8	
16	09/23/20 13:50:39.277.600	7f ff ff c8 8f ff ff c8	
17	09/23/20 13:50:39.281.900	9f ff ff c8 af ff ff c8	
18	09/23/20 13:50:39.289.900	bf ff ff c8 cf ff ff ff	
19	09/23/20 13:50:39.297.900	50 00 ff ff ff ff 10 1a	
20	09/23/20 13:50:39.350.800	f0 18 20 1e 73 10 30 1e	
21	09/23/20 13:50:39.359.200	00 14 40 13 41 e4 5f 7f	
22	09/23/20 13:50:39.367.500	ff c8 6f ff ff c8 7f ff	
23	09/23/20 13:50:39.375.500	ff c8 8f ff ff c8 9f ff	
24	09/23/20 13:50:39.390.900	ff c8 af ff ff c8 bf ff	
25	09/23/20 13:50:39.395.100	ff c8 cf 00 00 c7 50 00	
26	09/23/20 13:50:39.400.300	ff ff ff ff 10 1a fa c4	
27	09/23/20 13:50:39.452.200	20 1e 73 34 30 1d ff b0	
28	09/23/20 13:50:39.460.700	40 13 50 ac 5f ff ff c8	
29	09/23/20 13:50:39.475.600	6f ff ff c8 7f ff ff c8	
30	09/23/20 13:50:39.479.800	8f ff ff c8 9f ff ff c8	
31	09/23/20 13:50:39.486.100	af ff ff c8 bf ff ff c8	
32	09/23/20 13:50:39.494.500	cf 00 00 8f 50 00 ff ff	
33	09/23/20 13:50:39.502.800	ff ff	
34	09/23/20 13:50:39.546.600		30
35	09/23/20 13:50:39.548.500	10 1a fb 90 20 1e 73 34	
36	09/23/20 13:50:39.557.000	30 1e 00 8c 40 13 47 84	
37	09/23/20 13:50:39.573.400	5f 7f ff c8 6f ff ff c8	

Figure 4: Head Amp Data

## Data Structure

1. An example of a data capture is shown in file *Head Amp Capture.txt* supplied on the USB stick.
2. Note that there are 2 groups of data, DTE and DCE. This corresponds to data transmitted from and data received by the head amplifier respectively. The DCE data can be ignored as it serves no purpose.
3. The data transmitted by the head amp consists of 54 bytes sent every 100ms. The data is structured as detailed in Table 1 and Table 2.
4. Note that the terminator is 4 consecutive bytes of FF (hex). This can be used to align the data.
5. An extract from *Head Amp Capture.txt* is shown with the beam data highlighted in different colours. Beam 2 (Red), Beam 4 (Blue), Beam 1 (Green), Beam 3 (Orange)

Byte Number	Parameter	Collector	Beam Number	Full Scale (nA)	Value
0-3	Channel 1 Current	B	2	0.2	See Table 2
4-7	Channel 2 Current	D	4	20	
8-11	Channel 3 Current	A	1	20	
12-15	Channel 4 Current	C	3	0.1	
16-19	Channel 5 Current	E	5	Custom	
20-23	Channel 6 Current	F	6		
24-27	Channel 7 Current	G	7		
28-31	Channel 8 Current	H	8		
32-35	Channel 9 Current	Not connected			
36-39	Channel 10 Current	Not connected			
40-43	Channel 11 Current	J	9	Custom	
44-48	Channel 12 Current	K	10		
49-50	Gain Setting (Unused)	Not Applicable			Don't Care
50-53	Terminator	Not Applicable			FFFFFFFF(hex)

Table 1: Head Amp Data Structure

Byte Allocation (Big Endian Byte Order)	
D31..D28	D27..D0
Channel Number	Beam Current  Full Scale=DFFFC8(hex), 234,880,968(decimal) Overflow=FFFFFF(hex) D1..0 are always 0  Note: There is an electronic offset of 1% of full scale

Table 2: Byte Allocation

Line	MM/DD/YY	hh:mm:ss.ms.us	DTE(hex)	DCE(hex)
1	09/23/20	11:18:06.819.000	ff ff c8 bf ff ff c8 cf	
2	09/23/20	11:18:06.827.100	00 00 77 50 00 ff ff ff	
3	09/23/20	11:18:06.835.300	ff 10 1a ba 8c 20 1e 3e	
4	09/23/20	11:18:06.886.900	9c 30 1d b5 74 40 14 bf	
5	09/23/20	11:18:06.895.400	1c 5f 7f ff c8 6f ff ff	
6	09/23/20	11:18:06.908.700	c8 7f ff ff c8 8f ff ff	
7	09/23/20	11:18:06.912.900	c8 9f ff ff c8 af ff ff	
8	09/23/20	11:18:06.920.500	c8 bf ff ff c8 cf 00 00	
9	09/23/20	11:18:06.928.500	3f 50 00 ff ff ff ff 10	
10	09/23/20	11:18:06.981.600	1a a6 ec 20 1e 3f 20 30	
11	09/23/20	11:18:06.989.600	1d b5 d4 40 14 a1 2c 5f	
12	09/23/20	11:18:06.998.600	7f ff c8 6f ff ff c8 7f	
13	09/23/20	11:18:07.006.600	ff ff c8 8f ff ff c8 9f	
14	09/23/20	11:18:07.014.600	ff ff c8 af ff ff c8 bf	
15	09/23/20	11:18:07.027.800	ff ff c8 cf 00 00 07 50	
16	09/23/20	11:18:07.031.600	00 ff ff ff ff 10 1a af	
17	09/23/20	11:18:07.083.600	3c 20 1e 3e b8 30 1d b5	
18	09/23/20	11:18:07.091.600	7c 40 14 c7 34 5f 7f ff	
19	09/23/20	11:18:07.099.800	c8 6f ff ff c8 7f ff ff	
20	09/23/20	11:18:07.108.300	c8 8f ff ff c8 9f ff ff	
21	09/23/20	11:18:07.116.800	c8 af ff ff c8 bf ff ff	
22	09/23/20	11:18:07.125.300	c8 cf ff ff cf 50 00 ff	
23	09/23/20	11:18:07.133.700	ff ff ff	
24	09/23/20	11:18:07.145.700		0d
25	09/23/20	11:18:07.179.800	10 1a a6 10 20 1e 3e ec	
26	09/23/20	11:18:07.188.300	30 1d b6 78 40 14 ad c0	
27	09/23/20	11:18:07.205.300	5f ff ff c8 6f ff ff c8	
28	09/23/20	11:18:07.209.500	7f ff ff c8 8f ff ff c8	
29	09/23/20	11:18:07.213.800	9f ff ff c8 af ff ff c8	
30	09/23/20	11:18:07.221.800	bf ff ff c8 cf 00 00 97	
31	09/23/20	11:18:07.229.800	50 00 ff ff ff ff 10 1a	
32	09/23/20	11:18:07.283.000	a7 8c 20 1e 3e b0 30 1d	
33	09/23/20	11:18:07.291.400	b5 88 40 14 b0 80 5f 7f	
34	09/23/20	11:18:07.299.800	ff c8 6f ff ff c8 7f ff	
35	09/23/20	11:18:07.307.400	ff c8 8f ff ff c8 9f ff	
36	09/23/20	11:18:07.315.000	ff c8 af ff ff c8 bf ff	
37	09/23/20	11:18:07.328.800	ff c8 cf 00 00 5f 50 00	
38	09/23/20	11:18:07.332.800	ff ff ff ff 10 1a b2 1c	
39	09/23/20	11:18:07.384.100	20 1e 3d f8 30 1d b6 60	
40	09/23/20	11:18:07.392.700	40 14 c1 e4 5f 7f ff c8	

Figure 5: Example Data

NSF 23/09/20